



Mixing ratio	20-21 ml /100 gr
Mixing time by hand/Mixer	60" / 30"
Working time	9'
Final setting time (Vicat)	13'
Demoulding time	35' / 40'
Setting Expansion (2 h)	< 0.15%
Compression Strenght (1 h)	> 40 MPa

State of the art plaster. Thanks to its physical and mechanical characteristics, **PENTA ROCK** type V plaster, can be used to make models capable of duplicating prosthetics. It is very useful when a maximization of the master model is needed.

CODE	PACKAGING
DGP20 WT	Rectangular plastic drum kg 20
DGP5 WT	Conic plastic drum kg 5

Standard color: **WT**= Natural white

Color customization is available upon request.