



Mixing Ratio	22-23 ml /100 gr
Mixing Time by Hand / Mixer	60" / 30"
Working Time	7'
Final Setting Time (Vicat)	9'
Demoulding Time	35'
Setting Expansion (2 h)	< 0,10%
Compressive strength (1h)	> 35 MPa


Type IV plaster is easily workable with an excellent resistance to compression and abrasion.

MASTER ROCK is the ideal plaster for making denture models (mobile, total, partial, combined and skeletal) and master models.

CODE	PACKAGING
DGMS25 PK/GR/EW	White paper bag kg 25
DGM20 PK/GR/EW	Rectangular plastic drum kg 20
DGM25 PK/GR/EW	Conic plastic drum kg 25

Standard color: **PK**= Pink **GR**= Green **EW**= Extra white

Color customization is available upon request.



www.regeco.it
+39 0131 217398