



Mixing Ratio	26-27 ml /100 gr - 0.22 lbs
Mixing Time by Hand / Mixer	60" / 30"
Working Time	8'
Final Setting Time (Vicat)	10'
Demoulding Time	35'
Setting Expansion (2 h)	< 0.16%
Compressive strength (1h)	> 30 MPa

Type III plaster has a bright white color. It is ideal for orthodontic use.

KINDROCK offers smooth surfaces and easy polishing that will allow for the manufacturing of perfect “museum” models.

CODE	PACKAGING
DGKS25 EW	White paper bag kg 25
DGK20 EW	Rectangular plastic drum kg 20
DGK25 EW	Conic plastic drum kg 25

Standard color: **EW**= Extra white

Color customization is available upon request.